

Anti-Elmod2 Polyclonal Antibody For Wb, Elisa - Q8Iz81

Item Number: CM0030159

Introduction

ELMOD2 Rabbit Polyclonal Antibody for WB and ELISA. Detects human and mouse ELMOD2 (35kDa), a GTPase-activating protein involved in antiviral defense and Arf family GAP activity. Ideal for immunology and cell biology research.

[Learn More](#)

Parameter	Value
Product Name	ELMOD2 Rabbit pAb
Remarks/Alias	9830169G11Rik; ELMOD2
Species	Human
GeneID (Human)	255520
GeneID	255520
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of mouse ELMOD2 (Q8BGF6).)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8IZ81
Protein Weight	35kDa
Shipping	Ice bag