

Anti-ELN (Elastin) Rabbit Polyclonal Antibody For IHC-P, Elisa - P15502

Item Number: CM0017984

Introduction

CamelBio offers Anti-ELN (Elastin) Rabbit pAb, validated for IHC-P and ELISA. Targets human elastin (UniProt P15502), a 68kDa secreted extracellular matrix protein. Ideal for cardiovascular and connective tissue research.

[Learn More](#)

Parameter	Value
Product Name	ELN Rabbit pAb
Remarks/Alias	WS; WBS; SVAS; ADCL1; ELN
Species	Human
GeneID (Human)	2006
GeneID	2006
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-97 of human ELN (P15502).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	P15502
Protein Weight	68kDa
Shipping	Ice bag