

Anti-Elovl2 Polyclonal Antibody For Wb, Elisa - Q9Nxb9

Item Number: CM0018829

Introduction

Rabbit polyclonal antibody targeting human ELOVL2, a very long chain fatty acid elongase. Validated for WB and ELISA, cross-reacts with human and rat. Ideal for fatty acid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	ELOVL2 Rabbit pAb
Remarks/Alias	SSC2; ELOVL2
Species	Human
Gene ID (Human)	54898
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 247-296 of human ELOVL2 (NP_060240.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9NXB9
Protein Weight	35kDa
Shipping	Ice bag