

Anti-Elovl4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9GZr5

Item Number: CM0027697

Introduction

Rabbit polyclonal antibody targeting human ELOVL4, validated for Western blot, immunofluorescence/immunocytochemistry, and ELISA. Cross-reacts with mouse and rat. Ideal for studying very long chain fatty acid elongation and related diseases.

[Learn More](#)

Parameter	Value
Product Name	ELOVL4 Rabbit pAb
Remarks/Alias	ADMD; CT118; ISQMR; SCA34; STGD2; STGD3; ELOVL4
Species	Human
GeneID (Human)	6785
GeneID	6785
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ELOVL4 (NP_073563.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZR5
Protein Weight	37kDa
Shipping	Ice bag