

Anti-Elovl5 Polyclonal Antibody For Wb And Elisa - Q9Nyp7

Item Number: CM0023683

Introduction

Rabbit polyclonal antibody against human ELOVL5, validated for Western blot and ELISA. Targets very long chain fatty acid elongase 5, involved in lipid metabolism and linked to spinocerebellar ataxia.

[Learn More](#)

Parameter	Value
Product Name	ELOVL5 Rabbit pAb
Remarks/Alias	HELO1; SCA38; dj483K16.1; ELOVL5
Species	Human
GeneID (Human)	60481
GeneID	60481
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 247-326 of human ELOVL5 (NP_001229757.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NYP7
Protein Weight	35kDa
Shipping	Ice bag