

Anti-Elp2 Rabbit Polyclonal Antibody For Wb, Elisa - Q6la86

Item Number: CM0016562



Introduction

High-quality rabbit polyclonal antibody against human ELP2, validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for elongator complex and tRNA modification research.

[Learn More](#)

Parameter	Value
Product Name	ELP2 Rabbit pAb
Remarks/Alias	StIP; MRT58; SHINC-2; STATIP1; ELP2
Species	Human
Gene ID	55250
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-600 of human ELP2 (NP_001229808.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q6IA86
Protein Weight	93kDa
Shipping	Ice bag