

Anti-Elp5 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Te02

Item Number: CM0024464

Introduction

ELP5 Rabbit Polyclonal Antibody validated for WB and ELISA, reacts with human and mouse. Targets elongator complex protein 5 involved in tRNA modification and cell migration. Ideal for post-transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	ELP5 Rabbit pAb
Remarks/Alias	DERP6; MST071; HSPC002; MSTP071; C17orf81; ELP5
Species	Human
GeneID (human)	23587
GeneID	23587
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-316 of human ELP5 (NP_056177.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8TE02
Protein Weight	35kDa
Shipping	Ice bag