

Anti-Elspbp1 Polyclonal Antibody For Wb, If-P, Elisa - Q96Bh3

Item Number: CM0022096

Introduction

Rabbit polyclonal antibody against human ELSPBP1 (E12, HE12). Validated for WB, IF-P, ELISA; cross-reacts with mouse, rat. Detects secreted epididymal sperm-binding protein involved in fertilization.

[Learn More](#)

Parameter	Value
Product Name	ELSPBP1 Rabbit pAb
Remarks/Alias	E12; HE12; EL149; EDDM12; ELSPBP1
Species	Human
GeneID (Human)	64100
GeneID	64100
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 110-223 of human ELSPBP1 (NP_071425.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96BH3
Protein Weight	26kDa
Shipping	Ice bag