

Anti-Elspbp1 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Bh3

Item Number: CM0020156

Introduction

Rabbit polyclonal antibody raised against human ELSPBP1 protein. Validated for Western blot and ELISA applications. Aids research on sperm capacitation, fertilization, and reproductive biology.

[Learn More](#)

Parameter	Value
Product Name	ELSPBP1 Rabbit pAb
Remarks / Alias	E12, HE12, EL149, EDDM12, ELSPBP1
Species	Human
GeneID (Human)	64100
GeneID	64100
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-223 of human ELSPBP1 (NP_071425.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96BH3
Protein Weight	26kDa
Shipping	Ice bag