

Anti-Emb (Embigin/Gp70) Rabbit Polyclonal Antibody For Wb And Elisa - Q6Pcb8

Item Number: CM0012890



Introduction

Rabbit polyclonal antibody targeting human Embigin (EMB/GP70). Validated for WB and ELISA. Suitable for studies on motoneuron outgrowth, neuromuscular junctions, and monocarboxylate transporter localization.

[Learn More](#)

Parameter	Specification
Product Name	EMB Rabbit pAb
Remarks/Alias	GP70; EMB
Species	Human
GeneID (Human)	133418
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-260 of human EMB (NP_940851.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q6PCB8
Protein Weight	37kDa
Shipping	Ice bag