

Anti-Emc1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8N766

Item Number: CM0015244



Introduction

CamelBio supplies anti-EMC1 rabbit polyclonal antibody validated for WB, IHC-P, ELISA in human samples. Targets ER membrane protein complex subunit 1, ideal for IVD research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	EMC1 Rabbit pAb
Remarks/Alias	CAVIPMR; KIAA0090; EMC1
Species	Human
Gene ID (Human)	23065
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 885-965 of human EMC1 (NP_055862.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q8N766
Protein Weight	112kDa
Shipping	Ice bag