


Anti-Emc4 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q5J8M3

Item Number: CM0003830



Introduction

Rabbit monoclonal antibody against human ER membrane protein complex subunit 4 (EMC4, also PIG17/TMEM85). Validated for Western blot, IHC-P, and ELISA; cross-reacts with mouse and rat. For research use only.

[Learn More](#)

Parameter	Value
Product Name	EMC4 Rabbit mAb
Remarks/Alias	PIG17; TMEM85; EMC4
Species	Human
Gene ID (Human)	51234
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human EMC4 (NP_057538.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5J8M3
Protein Weight	15kDa/16kDa/20kDa
Shipping	Ice bag