


Anti-Emc8 Polyclonal Antibody For Wb, Elisa - O43402

Item Number: CM0013082



Introduction

EMC8 rabbit polyclonal antibody validated for Western blotting and ELISA, cross-reactive with human and mouse. Part of the EMC complex, essential for membrane protein insertion; suitable for ER and protein trafficking research.

[Learn More](#)

Parameter	Specification
Product Name	EMC8 Rabbit pAb
Remarks/Alias	NOC4; COX4NB; C16orf2; C16orf4; FAM158B; EMC8
Species	Human
Gene ID (Human)	10328
Gene ID	10328
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human EMC8 (NP_006058.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O43402
Protein Weight	24kDa
Shipping	Ice bag