

Anti-Eme1 Polyclonal Antibody For Wb And Elisa - Q96Ay2

Item Number: CM0030694



Introduction

Rabbit polyclonal antibody against EME1, suitable for Western blot and ELISA. Detects human EME1 (63 kDa) and cross-reacts with mouse, rat. Useful for studying DNA repair and genome stability.

[Learn More](#)

Parameter	Value
Product Name	EME1 Rabbit pAb
Remarks/Alias	MMS4L; SLX2A; EME1
Species	Human
GeneID (Human)	146956
GeneID	146956
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-320 of human EME1 (NP_689676.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96AY2
Protein Weight	63kDa
Shipping	Ice bag