


Anti-Emerin/Emd Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P50402

Item Number: CM0009963



Introduction

Rabbit monoclonal antibody against human Emerin/EMD, validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for nuclear envelope, chromatin, and muscular dystrophy research.

[Learn More](#)

Parameter	Value
Product Name	Emerin/EMD Rabbit mAb
Remarks/Alias	STA; EDMD; LEMD5; Emerin/EMD
Species	Human
Gene ID (Human)	2010
Gene ID	2010
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 150-254 of human Emerin/EMD (P50402))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50402
Protein Weight	29 kDa
Shipping	Ice bag