

Anti-Eml3 Polyclonal Antibody For Wb, Elisa - Q32P44

Item Number: CM0024380

Introduction

Rabbit polyclonal antibody targeting human EML3. Suitable for Western Blot and ELISA. Useful for studying mitotic spindle assembly, chromosome alignment, and microtubule nucleation.

[Learn More](#)

Parameter	Value
Product Name	EML3 Rabbit pAb
Remarks/Alias	ELP95; EMAP3; EMAP95; EML3
Species	Human
GeneID	256364
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-530 of human EML3 (NP_694997.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q32P44
Protein Weight	95kDa
Shipping	Ice bag