

Anti-Eml4 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Hc35

Item Number: CM0024528



Introduction

Rabbit polyclonal antibody against EML4, validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat EML4. Suitable for research on microtubule dynamics, mitosis, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	EML4 Rabbit pAb
Remarks/Alias	C2orf2; ELP120; EMAP-4; EMAPL4; ROPP120; EML4
Species	Human
GeneID (Human)	27436
GeneID	27436
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 360-560 of human EML4 (NP_061936.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HC35
Protein Weight	109kDa
Shipping	Ice bag