

Anti-Emr1 Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - Q14246

Item Number: CM0015967

Introduction

CamelBio's EMR1 Rabbit pAb (Q14246) is a polyclonal antibody for WB, IF/ICC, IF-P, ELISA in human, mouse, rat. Targets adhesion GPCR E1, involved in immune cell interactions and eosinophil function.

[Learn More](#)

Parameter	Value
Product Name	EMR1 Rabbit pAb
Remarks/Alias	EMR1; TM7LN3
Species	Human
GeneID (Human)	2015
GeneID	2015
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 477-596 of human EMR1 (NP_001965.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14246
Protein Weight	98kDa
Shipping	Ice bag