


Anti-Emr2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Uhx3

Item Number: CM0015567



Introduction

Rabbit polyclonal antibody targeting human EMR2 (ADGRE2, CD312), validated for WB and ELISA. Cross-reacts with mouse. Useful for studies on cell adhesion, chemotaxis, and mast cell degranulation.

[Learn More](#)

Parameter	Value
Product Name	EMR2 Rabbit pAb
Remarks/Alias	VBU; CD97; EMR2; CD312
Species	Human
Gene ID (Human)	30817
Gene ID	30817
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-500 of human EMR2 (NP_038475.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UHX3
Protein Weight	90kDa
Shipping	Ice bag