

## Anti-Emsy/C11Orf30 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q7Z589

Item Number: CM0030014



### Introduction

Rabbit polyclonal antibody against EMSY (C11orf30), a BRCA2-interacting transcriptional repressor involved in DNA repair and chromatin remodeling. Validated for WB, IF/ICC, ELISA in human and mouse.

[Learn More](#)

| Parameter        | Value                                                                                                                   |
|------------------|-------------------------------------------------------------------------------------------------------------------------|
| Product Name     | C11orf30 Rabbit pAb                                                                                                     |
| Remarks/Alias    | GL002; C11orf30                                                                                                         |
| Species          | Human                                                                                                                   |
| Gene ID          | 56946                                                                                                                   |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1153-1322 of human C11orf30 (NP_064578.2) |
| Source           | Rabbit                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                   |
| Application      | WB, IF/ICC, ELISA                                                                                                       |
| Cross Reactivity | Human, Mouse                                                                                                            |
| SWISS            | Q7Z589                                                                                                                  |
| Protein Weight   | 141kDa                                                                                                                  |
| Shipping         | Ice bag                                                                                                                 |