

Anti-Emx1 Polyclonal Antibody For Wb And Elisa - Q04741

Item Number: CM0024508

Introduction

Rabbit polyclonal antibody against human EMX1, suitable for WB and ELISA. Detects a 31 kDa protein in human cerebral cortex; cross-reacts with mouse. Supports developmental neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	EMX1 Rabbit pAb
Remarks/Alias	EMX1
Species	Human
GeneID	2016
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 70-190 of human EMX1 (NP_004088.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q04741
Protein Weight	31kDa
Shipping	Ice bag