

Anti-En2 Rabbit Polyclonal Antibody For Wb, Elisa - P19622

Item Number: CM0022891

Introduction

Rabbit polyclonal antibody against human EN2, validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for studies on homeobox protein engrailed-2 in development and autism spectrum disorder research.

[Learn More](#)

Parameter	Value
Product Name	EN2 Rabbit pAb
Remarks/Alias	EN2
Species	Human
GeneID (Human)	2020
GeneID	2020
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 250-333 of human EN2 (NP_001418.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19622
Protein Weight	34kDa
Shipping	Ice bag