

Anti-Enah Polyclonal Antibody For Wb, If/Icc, Elisa - Q8N8S7

Item Number: CM0021336



Introduction

ENAH Rabbit pAb is a polyclonal antibody validated for WB, IF/ICC, and ELISA, targeting human ENAH (also known as MENA). Cross-reacts with mouse and rat. Optimal for cytoskeleton remodeling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	ENAH Rabbit pAb
Remarks/Alias	ENA, MENA, NDPP1, ENAH
Species	Human
GeneID (Human)	55740
GeneID	55740
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 412-591 of human ENAH (NP_001008493.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N8S7
Protein Weight	67kDa
Shipping	Ice bag