

Anti-Endo G Rabbit Monoclonal Antibody For Wb, Elisa - Q14249

Item Number: CM0008855

Introduction

Recombinant rabbit monoclonal antibody targeting human Endonuclease G, mitochondrial (ENDO G). Validated for Western blot and ELISA. Detects human, mouse, and rat Endo G. Suitable for apoptosis and mitochondrial DNA research.

[Learn More](#)

Parameter	Value
Product Name	Endo G Rabbit mAb
Remarks/Alias	ENDO G; endonuclease G; Endo G
Species	Human
Gene ID	2021 (human)
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 198-297 of human Endo G (Q14249)
Source	Rabbit
Category	Monoclonal Antibody
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14249
Protein Weight	33kDa
Shipping	Ice bag