

Anti-Endogl1/Engl Rabbit Monoclonal Antibody For Wb And Elisa - Q9Y2C4

Item Number: CM0002780



Introduction

Rabbit monoclonal antibody targeting human ENDOGL1/ENGL (EXOG), a mitochondrial nuclease. Ideal for Western blot (WB) and ELISA applications. Suitable for investigating mitochondrial DNA maintenance and apoptosis. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	ENDOGL1/ENGL Rabbit mAb
Remarks/Alias	ENGL; ENGLA; ENGLB; ENGL-a; ENGL-b; ENDOGL1; ENDOGL2; ENDOGL1 / ENGL
Species	Human
Gene ID (Human)	9941
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 269-368 of human ENDOGL1 / ENGL (Q9Y2C4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y2C4
Protein Weight	41kDa
Shipping	Ice bag