

Anti-Endomucin Polyclonal Antibody For Wb, Elisa - Q9ULc0

Item Number: CM0028016

Introduction

Endomucin rabbit polyclonal antibody validated for Western blot and ELISA. Reacts with human and mouse endomucin. Ideal for studying endothelial cell function, vascular biology, and extracellular matrix interactions.

[Learn More](#)

Parameter	Value
Product Name	Endomucin Rabbit pAb
Remarks/Alias	EMCN2; MUC14; Endomucin
Species	Human
Gene ID (Human)	51705
Gene ID	51705
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 19-135 of human Endomucin (NP_057326.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9ULC0
Protein Weight	27kDa
Shipping	Ice bag