

Anti-Endonuclease G (Endo G) Rabbit Polyclonal Antibody For Wb, Elisa - Q14249

Item Number: CM0016515

Introduction

Rabbit polyclonal antibody against human Endonuclease G (Endo G) validated for WB and ELISA; cross-reacts with mouse and rat. Useful for studying mitochondrial endonuclease function in DNA repair, autophagy, and stress responses.

[Learn More](#)

Parameter	Value
Product Name	Endo G Rabbit pAb
Remarks/Alias	ENDOG; Endo G
Species	Human
GeneID (human)	2021
GeneID	2021
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 49-297 of human Endo G (NP_004426.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14249
Protein Weight	33kDa
Shipping	Ice bag