


Anti-Endothelin 1 (Et1, Ppet1) Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P05305

Item Number: CM0015619



Introduction

CamelBio offers rabbit anti-Endothelin 1 polyclonal antibody (ET1, PPET1) validated for WB, IF/ICC, IHC-P, and ELISA. Reacts with human, mouse, and rat samples. Ideal for cardiovascular and placental research.

[Learn More](#)

Parameter	Value
Product Name	Endothelin 1 Rabbit pAb
Remarks/Alias	ET1; QME; PPET1; ARCND3; HDLCQ7; Endothelin 1
Species	Human
GeneID (human)	1906
GeneID	1906
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-212 of human Endothelin 1 (NP_001946.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05305
Protein Weight	24kDa
Shipping	Ice bag