

Anti-Engrailed 1 Monoclonal Antibody For Wb - Q05925

Item Number: CM0002939

Introduction

Anti-Engrailed 1 rabbit monoclonal antibody validated for WB and ELISA applications. Recognizes human, mouse, and rat Engrailed 1. Suitable for studies on limb development and dorsal-ventral patterning.

[Learn More](#)

Parameter	Value
Product Name	Engrailed 1 Rabbit mAb
Remarks/Alias	ENDOVESLB; Engrailed 1
Species	Human
GeneID (human)	2019
GeneID	2019
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 293-392 of human Engrailed 1 (Q05925).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q05925
Protein Weight	40kDa
Shipping	Ice bag