


Anti-Eno1 Rabbit Monoclonal Antibody For Wb, Elisa - P06733

Item Number: CM0002911



Introduction

ENO1 Rabbit mAb (P06733) detects alpha-enolase in human, mouse, and rat samples. Validated for Western blot and ELISA, targeting a 47 kDa protein. Suitable for glycolysis, cancer, and fibrinolytic studies.

[Learn More](#)

Parameter	Value
Product Name	ENO1 Rabbit mAb
Remarks/Alias	NNE; PPH; MPB1; ENO1L1; ENO1-IT1; HEL-S-17; ENO1
Species	Human
Gene ID (Human)	2023
Gene ID	2023
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ENO1 (P06733)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06733
Protein Weight	47kDa
Shipping	Ice bag