

Anti-Eno1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P06733

Item Number: CM0022133



Introduction

Validate ENO1 expression in human, mouse, and rat samples with this rabbit polyclonal antibody. Suitable for western blot, immunofluorescence, immunoprecipitation, and ELISA; recombinant immunogen ensures specificity for glycolysis and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	ENO1 Rabbit pAb
Remarks/Alias	NNE; PPH; MPB1; ENO1L1; ENO1-IT1; HEL-S-17; ENO1
Species	Human
Gene ID (Human)	2023
Gene ID	2023
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-434 of human ENO1 (NP_001419.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06733
Protein Weight	47kDa
Shipping	Ice bag