

Anti-Eno2 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P09104

Item Number: CM0017488



Introduction

Rabbit polyclonal antibody against human ENO2 (NSE), validated for WB, IF/ICC, IP, ELISA. Cross-reacts with mouse and rat. Ideal for neuroscience and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	ENO2 Rabbit pAb
Remarks/Alias	NSE; HEL-S-279; ENO2
Species	Human
GeneID (Human)	2026
GeneID	2026
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 195-434 of human ENO2 (NP_001966.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09104
Protein Weight	47kDa
Shipping	Ice bag