


Anti-Eno3 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P13929

Item Number: CM0004482



Introduction

Rabbit monoclonal antibody against human ENO3 (Beta-enolase, MSE), validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for muscle metabolism and glycogen storage disease research.

[Learn More](#)

Parameter	Value
Product Name	ENO3 Rabbit mAb
Remarks/Alias	MSE; GSD13
Species	Human
GeneID (Human)	2027
GeneID	2027
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human ENO3 (NP_001967.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13929
Protein Weight	42kDa/44kDa/47kDa
Shipping	Ice bag