


Anti-Eno3 Rabbit Polyclonal Antibody For Wb And Elisa - P13929

Item Number: CM0031014



Introduction

Rabbit polyclonal antibody against human ENO3, validated for WB and ELISA. Cross-reacts with human, mouse, and rat beta-enolase. Ideal for glycolysis and muscle development research.

[Learn More](#)

Parameter	Value
Product Name	ENO3 Rabbit pAb
Remarks/Alias	MSE; GSD13; ENO3
Species	Human
GeneID (Human)	2027
GeneID	2027
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human ENO3 (NP_001967.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13929
Protein Weight	47kDa
Shipping	Ice bag