


Anti-Eno3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P13929

Item Number: CM0020225



Introduction

Rabbit polyclonal antibody against human ENO3 (beta-enolase, muscle-specific enolase). Validated for WB, IF/ICC, ELISA in human, mouse, rat samples. Ideal for studying glycolysis, muscle development, and glycogen storage disease.

[Learn More](#)

Parameter	Value
Product Name	ENO3 Rabbit pAb
Remarks/Alias	MSE; GSD13; ENO3
Species	Human
GeneID (Human)	2027
GeneID	2027
Immunogen	Recombinant Protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human ENO3 (NP_001967.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13929
Protein Weight	47kDa
Shipping	Ice bag