

Anti-Enos Monoclonal Antibody For Wb, Ihc-P, Elisa - P29474

Item Number: CM0003456



Introduction

Rabbit monoclonal antibody specific for human eNOS (NOS3), validated for WB, IHC-P, and ELISA applications. Cross-reacts with mouse and rat. Targets endothelial nitric oxide synthase, a key enzyme in vascular biology.

[Learn More](#)

Parameter	Value
Product Name	eNOS Rabbit mAb
Remarks/Alias	eNOS; ECNOS
Species	Human
GeneID (human)	4846
GeneID	4846
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1035-1110 of human eNOS (NP_000594.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29474
Protein Weight	133kDa
Shipping	Ice bag