

Anti-Enox2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q16206

Item Number: CM0026323

Introduction

Rabbit polyclonal ENOX2 antibody for WB, IF/ICC, ELISA. Reacts with human, mouse, rat ENOX2 (APK1, tNOX), a tumor-associated cell membrane oxidase. Ideal for cancer research.

[Learn More](#)

Parameter	Value
Product Name	ENOX2 Rabbit pAb
Remarks/Alias	APK1, tNOX, COVA1, ENOX2
Species	Human
Gene ID (Human)	10495
Gene ID	10495
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 411-610 of human ENOX2 (NP_872114.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16206
Protein Weight	70kDa
Shipping	Ice bag