

Anti-Enpp1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - G3X9S2

Item Number: CM0021454

Introduction

Rabbit polyclonal antibody against mouse Enpp1 (Ectonucleotide pyrophosphatase/phosphodiesterase family member 1), validated for WB, IF/ICC, and ELISA. Cross-reacts with human and rat. UniProt accession G3X9S2.

[Learn More](#)

Parameter	Value
Product Name	Enpp1 Rabbit pAb
Remarks/Alias	Pca, ttw, twy, M6S1, NPP1, Npps, PC-1, Ly-41, Pca-1, Pdnp1, CD203c, E-NPP1, E-NPP 1, 4833416E15Rik, Enpp1
Species	Mouse
GeneID (Human)	5167
GeneID (Mouse)	18605
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 580-730 of mouse Enpp1 (NP_032839.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	G3X9S2
Protein Weight	103kDa
Shipping	Ice bag