


Anti-Enpp3 Polyclonal Antibody For Wb, Elisa - 014638

Item Number: CM0018525



Introduction

ENPP3 rabbit polyclonal antibody validated for WB and ELISA. Targets ectonucleotide pyrophosphatase/phosphodiesterase 3 (CD203c), involved in nucleotide metabolism and immune regulation. Cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	ENPP3 Rabbit pAb
Remarks/Alias	B10; NPP3; PDNP3; CD203c; PD-IBETA; ENPP3
Species	Human
GeneID (Human)	5169
GeneID	5169
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 646-875 of human ENPP3 (NP_005012.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14638
Protein Weight	100kDa
Shipping	Ice bag