


Anti-Enpp3/Cd203C Rabbit Monoclonal Antibody For Wb, If/Icc, Fc, Elisa - O14638

Item Number: CM0007352



Introduction

ENPP3/CD203c rabbit monoclonal antibody validated for WB, IF/ICC, FC, ELISA. Targets human ectonucleotide pyrophosphatase/phosphodiesterase 3, a key regulator of extracellular nucleotide metabolism and innate immunity. Ideal for immunology and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	ENPP3/CD203c Rabbit mAb
Remarks/Alias	B10; NPP3; PDNP3; CD203c; PD-IBETA; ENPP3/CD203c
Species	Human
Gene ID (Human)	5169
Gene ID	5169
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 551-875 of human ENPP3/CD203c (NP_005012.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	O14638
Protein Weight	100kDa
Shipping	Ice bag