

Anti-Enpp5 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Uja9

Item Number: CM0009477



Introduction

Rabbit monoclonal antibody targeting human ENPP5, suitable for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Supports research in neuronal communication and ectonucleotidase function.

[Learn More](#)

Parameter	Value
Product Name	ENPP5 Rabbit mAb
Remarks/Alias	NPP5; NPP-5; ENPP5
Species	Human
GeneID (human)	59084
GeneID	59084
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 378-477 of human ENPP5 (Q9UJA9).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UJA9
Protein Weight	55kDa
Shipping	Ice bag