

Anti-Ensa/Arpp-19E Polyclonal Antibody For Wb And Elisa - O43768

Item Number: CM0017266

Introduction

Rabbit polyclonal antibody against human alpha-endosulfine (ENSA/ARPP-19e), suitable for WB and ELISA, detects mouse protein. Ideal for studying PP2A inhibition during mitosis and insulin secretion pathways.

[Learn More](#)

Parameter	Value
Product Name	ENSA Rabbit pAb
Remarks/Alias	ARPP-19e; ENSA
Species	Human
GeneID (Human)	2029
GeneID	2029
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ENSA (NP_996929.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O43768
Protein Weight	13kDa
Shipping	Ice bag