

Anti-Ent2/Slc29A2 Rabbit Monoclonal Antibody For Ihc-P, If/Icc, Elisa - Q14542

Item Number: CM0002966

Introduction

Rabbit monoclonal antibody against human ENT2/SLC29A2 (Q14542), suitable for IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying nucleoside transport and equilibrative nucleoside transporter function.

[Learn More](#)

Parameter	Value
Product Name	ENT2/SLC29A2 Rabbit mAb
Remarks/Alias	ENT2; DER12; HNP36; ENT2/SLC29A2
Species	Human
GeneID (Human)	3177
GeneID	3177
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ENT2/SLC29A2 (Q14542).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14542
Protein Weight	50kDa
Shipping	Ice bag