

Anti-Entpd1/Cd39 Polyclonal Antibody For Wb, If, Ihc-P - P49961

Item Number: CM0026064

Introduction

ENTPD1/CD39 Rabbit pAb targets ectonucleoside triphosphate diphosphohydrolase 1, a key regulator of purinergic signaling. Validated for WB, IF/ICC, IF-P, IHC-P, ELISA with cross-reactivity to human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	ENTPD1/CD39 Rabbit pAb
Remarks/Alias	CD39; SPG64; ATPDase; NTPDase-1; ENTPD1/CD39
Species	Human
GeneID (Human)	953
GeneID	953
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human ENTPD1/CD39 (NP_001157655.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49961
Protein Weight	58kDa
Shipping	Ice bag