

Anti-Entpd5 Monoclonal Antibody For Wb, Ihc-P, Elisa - O75356

Item Number: CM0002147

Introduction

Rabbit monoclonal antibody against ENTPD5, validated for WB, IHC-P, ELISA in human, mouse, rat. Targets a 48 kDa ER nucleoside diphosphatase involved in glycoprotein quality control.

[Learn More](#)

Parameter	Specification
Product Name	ENTPD5 Rabbit mAb
Remarks/Alias	PCPH; CD39L4; NTPDase-5; ENTPD5
Species	Human
GeneID (Human)	957
GeneID	957
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human ENTPD5 (O75356).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75356
Protein Weight	48kDa
Shipping	Ice bag