

Anti-Entpd5 Rabbit Monoclonal Antibody For Wb, Elisa - 075356

Item Number: CM0002478

Introduction

Rabbit monoclonal antibody against human ENTPD5 (CD39L4/PCPH), suitable for Western blot and ELISA, with cross-reactivity to mouse and rat. Ideal for studying ER UDPase function in glycoprotein folding and quality control.

[Learn More](#)

Parameter	Value
Product Name	ENTPD5 Rabbit mAb
Remarks/Alias	PCPH; CD39L4; NTPDase-5; ENTPD5
Species	Human
GeneID	957
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 300-400 of human ENTPD5 (075356).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	075356
Protein Weight	48kDa
Shipping	Ice bag