

Anti-Entpd6 Polyclonal Antibody For Wb, Ihc-P, Elisa - 075354

Item Number: CM0023026

Introduction

ENTPD6 Rabbit pAb is a polyclonal antibody raised against the human ENTPD6 protein (aa 61-300). Validated for WB, IHC-P, and ELISA, it recognizes human, mouse, and rat ENTPD6, a Golgi membrane nucleoside diphosphatase involved in nucleotide hydrolysis and glycosylation.

[Learn More](#)

Parameter	Value
Product Name	ENTPD6 Rabbit pAb
Remarks/Alias	CD39L2; IL6ST2; IL-6SAG; NTPDase-6; dj738P15.3; ENTPD6
Species	Human
GeneID (Human)	955
GeneID	955
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 61-300 of human ENTPD6 (NP_001304870.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075354
Protein Weight	53kDa
Shipping	Ice bag