

Anti-Eomes Polyclonal Antibody For Wb, Elisa - 095936

Item Number: CM0018759

Introduction

Rabbit polyclonal antibody targeting human EOMES (TBR2). Validated for WB and ELISA; reacts with human, mouse, and rat. Ideal for studies in developmental biology, neurobiology, and immune response regulation.

[Learn More](#)

Parameter	Value
Product Name	EOMES Rabbit pAb
Remarks/Alias	TBR2; EOMES
Species	Human
GeneID (human)	8320
GeneID	8320
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 463-686 of human EOMES (NP_005433.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	095936
Protein Weight	73kDa
Shipping	Ice bag