

Anti-Ep300 Polyclonal Antibody For Ihc-P, Elisa - Q09472

Item Number: CM0020723



Introduction

Rabbit polyclonal antibody against human histone acetyltransferase p300 (EP300). Validated for IHC-P and ELISA, cross-reacts with mouse and rat. Supports chromatin biology and transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	EP300 Rabbit pAb
Remarks/Alias	p300; KAT3B; MKHK2; RSTS2; EP300
Species	Human
Gene ID (Human)	2033
Gene ID	2033
Immunogen	Recombinant protein corresponding to a sequence within amino acids 940-1030 of human EP300 (NP_001420.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q09472
Protein Weight	264kDa
Shipping	Ice bag