

Anti-Ep300 Polyclonal Antibody For Ihc-P, If/Icc, Elisa - Q09472

Item Number: CM0017427



Introduction

Rabbit polyclonal EP300 antibody validated for IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat EP300. Suitable for chromatin remodeling and transcriptional coactivation studies.

[Learn More](#)

Parameter	Value
Product Name	EP300 Rabbit pAb
Remarks/Alias	p300; KAT3B; MKHK2; RSTS2; EP300
Species	Human
Gene ID (Human)	2033
Gene ID	2033
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human EP300 (NP_001420.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q09472
Protein Weight	264kDa
Shipping	Ice bag